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Patient 3D Dosimetry

Measurement Evaluation **Dose Evaluation** Side by side 3 plane dose Statistics

Absolute dose Relative reference dose maximum Relative entered value 6977.7 [cGy]

Select reference dose Tps 01_HnN66vmat. TPS dose

Select evaluation dose 01_HnN66vmat. Reconstructed with HM_0, Machine: Clinac4_Tallin_2 NON-CLINICAL

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01_HnN66vmat. TPS dose 01_HnN66vmat. Reconstructed with HM_0, Machine: Clina...

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2D Plan Verification


Patient 3D Dosimetry

[Not connected] Administrator, Administrator


17:19 07.05.2015

Evaluación 3D DVH solución 2-1


TPS

Compass
Computed (CC)

Compass
Reconstructed (CR)


↔
 Δ, γ

↔
 Δ, γ

98% of beam information
missed. Point measurements
limit predictability*

Dolphin
MatriXX Ion Chamber

Accurate and high-sensitivity measurements

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


myQA®


LWS Beck VMAT.0: TPS dose

LWS Beck VMAT.0: Reconstructed with 18.01.2017 16:51:09, Machine: Versa26jul16

Imaging system: AS Open20, Outline: HC_Mensch, Dose grid resolution [cm]: [0,4 0,4 0,4]

Instant check on the actual delivered fluence

The screenshot displays the COMPASS software interface, which is used for monitoring and evaluating the actual delivered fluence. The interface consists of several overlapping windows, each showing a detector response histogram. These histograms plot 'Amount [%]' on the y-axis against 'Difference [%] of' on the x-axis. The histograms show the distribution of response differences for various measurement intervals and beam coordinates.

Key parameters visible in the histograms include:

- Beam: ARC 1
- Gantry: 37.23°
- Measurement interval: 17.09
- Beam coordinate: 18.24-20.5 cm
- Beam coordinate: 18.33-12.47 cm
- Beam coordinate: 14.68-0.80 cm
- Beam coordinate: 24.55-9.64 cm

A central dialog box titled "QuickCheck evaluation completed" is overlaid on the interface. It features a green checkmark icon and the text "Overall evaluation status is Passed." Below this, it asks "Would you like to show detail results immediately?" with "Yes" and "No" buttons. The dialog also indicates "Time to close: 2 s".

In the background, a larger histogram titled "Selected Pixel Count per Degree in S" is visible, plotting "Counts/Degree" on the y-axis against "Degree" on the x-axis. This histogram compares "Predicted" (blue line) and "Measured" (red line) data, showing a significant peak around 60 degrees.

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File Device Computation Tools Help

Patient 3D Dosimetry

Measurement Evaluation Dose Evaluation Side by side 3 plane dose Statistics

Absolute dose Relative reference dose maximum Relative entered value 3010.5 [cGy]

Select reference dose Tps COTE G.O: TPS dose

Select evaluation dose COTE G.O: Reconstructed with 26.08.2016 15:47:08, Machine: Agility_Rennes_1

Position: 17,87 1,00 -45,20 cm
CT: - HU
Density: -
Dose: -

Position: 17,87 1,00 -45,20 cm
CT: - HU
Density: -
Dose: -

Position: 17,87 1,00 -45,20 cm
CT: - HU
Density: -
Gamma: -

Beam: 10A0 (Gantry: 340,00° Collimator: 0,00° Couch: 0,00°)
Transversal: 1,00 cm
Slice 68/192

Beam: 10A0 (Gantry: 340,00° Collimator: 0,00° Couch: 0,00°)
Transversal: 1,00 cm
Slice 68/192

Beam: 10A0 (Gantry: 340,00° Collimator: 0,00° Couch: 0,00°)
Transversal: 1,00 cm
Slice 68/192

Regions of interest

ROi	Matl	Type
<input checked="" type="checkbox"/>	p drt	Organ
<input checked="" type="checkbox"/>	p pcl	Organ
<input checked="" type="checkbox"/>	moel	Organ
<input checked="" type="checkbox"/>	GTV	Organ
<input checked="" type="checkbox"/>	Régl	External
<input checked="" type="checkbox"/>	ITV	Organ
<input checked="" type="checkbox"/>	PTV	Organ
<input checked="" type="checkbox"/>	Coel	Organ
<input checked="" type="checkbox"/>	estor	Organ
<input checked="" type="checkbox"/>	Cext	Organ
<input checked="" type="checkbox"/>	table	Support
<input checked="" type="checkbox"/>	table	Support
<input checked="" type="checkbox"/>	Can	Organ
<input checked="" type="checkbox"/>	CON	Organ

Statistics

Dvh Default Statistics Clinical goals Comparative Goals

COTE G.O: TPS dose COTE G.O: Reconstructed with 26.08.2016 15:47:08, Mach...

DVH

Distance criteria: 0,3 cm Dose percentage: 3%

Position: 17,87 1,00 -45,20 cm
CT: - HU
Density: -
Gamma: -

Beam: 10A0 (Gantry: 340,00° Collimator: 0,00° Couch: 0,00°)
Transversal: 1,00 cm
Slice 68/192

Partnered with Iba RaySearch Laboratories


Detector Setup 2D Plan Verification Patient 3D Dosimetry

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Godart et al. Med. Phys. Biol, 56 (2011)

Implementation of Stereotactic Ablative Body Radiotherapy (SABR) for early stage medically inoperable lung cancer: An experience from a medium size clinical centre, UKRO conference 2015

Stereotactic Radiation Therapy (SRS/SBRT) pre-treatment QA: two different approaches, 3rd ESTRO Forum (Barcelona 2015).

Lung SBRT treatment verification with a transmission detector system, ESTRO2017

■ ¿qué es myQA?

Agenda_DS... Dropbox SaralC82

BistroMath MagicMaX OmniPro Advance 1.2 SARRP presentatio... expenses

DigiPhant Deskshare Bulk Rename Utility Lynx 2D QPimRT1.7 Siemens Artist Loca... PT Brochures.rar

HyperCam 2 - RadioModels Lynx2DIm... SCH_Sales (SCHFILE1)... RE COMPA... Agenda_DS... Template DIN protoc...

OmniProIn... OmniPro... 7.5 Microsoft Outlook 2010 ZETTERFAS... Clinical implement... folder.zip Template DIN protoc...

wetransfer... Treviso DOPA - Shortcut Clinical Analysis Kagoshima... DIN protocol for myQA...

Galway COMPASS AAPM_Hig... Siemens_O... PMB Copy of CCU 10FFF_test...

AAPM characteriz... transmissio... COMPASS, KFMC wor... RT-BR-E-D... Copy of P1 Dosi Packa... Galaway.rs...

Poznan_1_... Neso_Ennio... PT Brochures Workshop Final IBA.pdf DOSI_PT_P... Galaway_s...

RT-BR-E-M... H+N LET protons.pdf HIN OTP.20151... P1 Dosi Package_M... Treviso Agil_S.rsbak

folder myQA versi... MRI linac.pdf Copy of CCU upgr. CUF... Service_Trai...

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









myQA Modules


Generic Tests




Check lists, pass/fail, functional tests, numeric tests...



Dosimetric Tests

MLC/VMAT tests

EPID/KV/CBCT Tests/WL test

myQA Plug-ins

myQA Platform

- myQA Platform
- myQA Patients
- myQA Machines
- myQA FastTrack
- myQA Accept

Equipment Setup

User Administration

Options

Data Import

Compare

Welcome Poznan

You are assigned to Poznan Oncology Center

myQA RSS Feed

[Register now for myQA User Meeting on November 23+ 24 in Germany](#) 13.09.2017
 Register now for a myQA User Meeting on November 23 and 24 in our ICC in Schwarzenbruck, Germany. Want to learn more about this event? Then just click above link.

[myQA Patients - Technical Notes published on Support Portal](#) 05.09.2017
 Information for users of myQA Patients: two relevant technical notes have been published on our Support Portal. These are: "myQA Patients_TN008_170808 01_myQA Patients_Normalization lost" and "myQA Patients_TN009_170818 01_Align Cube on Isocenter". The latter applies to OP-ImRT+ as well. The technical notes are available from the "Useful Technical Materials" and connected to the product "myQA Patients".

[Release of myQA 2017_001 - Distribution via Support Portal](#) 22.06.2017
 We are happy to announce the release of myQA® Version 2017-001 (Build 2.8.18). This version will be distributed via the IBA Support Portal. Eligible customers who have registered themselves in the portal will receive a download ticket in the coming days by email from mailer@iba-portal.com. Key features of this release include: Support of Windows 10 and Virtual Machines, simplified initial setup and significantly speed up "qoing"

Training

TRAINING COURSES FOR HEALTHCARE PROFESSIONALS

ICC [INTERNATIONAL COMPETENCE CENTER]

Service

SUPPORT PORTAL

Your 24/7 access to IBA service

CAREprogram

myQA Platform

Main

Start Page

Equipment Setup

User Administration

Options

Data Import

Compare

Welcome Pozna

- ▶ Demo Clinic_1
 - ▶ Wielkopolskie Centrum Onkologii
 - ▶ Diagnostic Department_Poznan
 - ▶ CT Room
 - ▶ CT Poznan
 - ▶ Ortodontic
 - ▶ Ortodontic_PZ
 - ▶ RT Department_Poznan
 - ▶ Cyber Knife Room
 - ▶ Cy PZ 01
 - ▶ GR Room
 - ▶ GR_Lynx
 - ▶ GR1_P1
 - ▶ Linac Room 1
 - ▶ TB01
 - ▶ Linac Room 2
 - ▶ Linac Room 2 ST
 - ▶ TB02
 - ▶ Room 3
 - ▶ R3 Halcyon
 - ▶ Tomo Room
 - ▶ Tomo
 - ▶ TPS room
 - ▶ RS

myQA RSS Feed

[Register now for myQA User M](#)
 Register now for a myQA User learn more about this event? TI

[myQA Patients - Technical Not](#)
 Information for users of myQA These are: "myQA Patients_TNC Patients_TN009_170818 01_Aliq are available from the "Useful 1

[Release of myQA 2017 001 - D](#)
 We are happy to announce the distributed via the IBA Support receive a download ticket in th include: Support of Windows 10 and virtual machines, simplified initial setup and significantly speed up...

Edit

Wielkopolskie Centrum Onkologii

Details

Clinic code:

Address: Ul. Gabary 15
61-866 Poznan, Poland

Operator:

Phone:

Description:

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myQA - 2017-002 - Not for clinical use!

myQA Machines

Main | Test Run

Clinics | Add Unscheduled | Skip | Finish | Skip | Reset | Override Status | Do not group | Group by Category | Group by Status | Refresh

Locations | Active Tasks | Task | Test | View | Data

Iba

Poznan

Machines

CT Poznan	<none>	(0)
Cy PZ 01	<none>	(0)
TB01	<none>	(0)
TB02	Poznan QA Daily	(2)
Tomo	<none>	(0)

Active Tasks

I - Daily QA	tomorrow	(0/13)
II - special Tests...	tomorrow	(0/11)

I - Daily QA due tomorrow

013a_MLC accuracy	MLC	Start
013b_Picket Fence test Dynamic	MLC	Start
013c_Picket Fence Dinamic with Error	MLC	Start
014a_VMAT leaves speed	MLC	Start
014b_Gantry speed	MLC	Start
015a_CBCT Contrast check	CBCT kv/MV	Start
015b_CBCT HU constancy	Other Energy: n/a <input checked="" type="checkbox"/> Share Results HU Constancy test was cancelled.	Start
016a_EPID/KV scaling	Planar MV (EPID)	Start
016b_EPID/KV spatial resolution	Planar kv/MV	Start
016c_EPID/KV Contrast	Planar kv/MV	Start
016D1_Las Vegas Visual Test	Planar MV (EPID)	Pass Warning Fail
016D2_Las Vegas Visual Test	Planar MV (EPID)	Pass Warning Fail

Task

Description

From AAPM TG-142:
"The daily tests include parameters that can affect dose to the patient by dosimetric (output constancy) or geometric (lasers, optical distance indicator, field size) means. The daily safety tests still include audiovisual monitoring of the patient and testing of the door interlock. With respect to EPID and kV imaging, the operation and functionality are tested daily, as well as collision interlocks".

The tests are the ones listed in the Table I, Table IV and Table VI of the AAPM TG-142 Report.

Note Edit

Due date changed from 12.10.2017 to 06.10.2017 by Poznan

Test: 1.01 Output constancy wi...

Description

Poznan QA Procedure for Daily output. use the plate with he Farmer chamber. Register the dose directly in myQA or on the xls sheet

Attachments


[IMG_2131\[1\].JPG](#)
[IMG_2136\[1\].JPG](#)
[myQA-white-paper_Tg-142_Daily-Generac-QA_R1_2015-04-20](#)

Note Edit

<none>

Reference Date

Enter date 05

myQA: TG-66

myQA - 2016-002 - Not for clinical use!

Administrator

myQA Machines

Test Setup

Test Run

Test Repository

Locations Edit Clipboard Item Protocol Template Protocol Data View

Machines & Templates

- Protocol Templates
 - CT-dedicated <none>
 - Elekta Demo Protocol
 - P1 <none>

Protocols: CT-dedicated

- TG-66
 - Annual QA
 - 1.0 Gantry Tilt Accuracy (...)
 - 1.1 Table Indexing and Po...
 - 1.2 Scan Localization (Me...
 - 2.0 CTDI (Dosimetry)
 - Daily QA
 - Monthly QA

Test Definition

Name: 1.0 Gantry Tilt Accuracy

Description: This test will verify the accuracy of the gantry tilt indicators

1. Setup gantry with digital level or inclinometer on face of gantry near control module
2. Tilt the gantry to the desired angle and record the reading
3. Bring the gantry back to 0 degrees and record the reading

If you lasers are mounted directly in front of the gantry that would inhibit gantry tilt do not proceed with this test
 If you have a ceiling mounted power injector make sure it is clear of the gantry prior to start of test

Acceptance Criteria

Name	Expected / Ref	Warn	Fail	Type	Relative
10 degree Offset	10,0 °	± 0,5 °	± 1,0 °	Tolerance	No
20 degree offset	20,0 °	± 0,5 °	± 1,0 °	Tolerance	No
Return to 0 degree	0,0 °	± 0,5 °	± 1,0 °	Tolerance	No

Mechanical
Numeric

Stop Watch

17:04
29.08.2016

IBA Dosimetry - myQA Cockpit

http://localhost/#/

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- Agility 4578
- I - Daily QA
- II - Weekly QA
- III - Monthly QA
- IV - Annual QA
- Artisle 123
- Task CBCT
- Task Dosi
- Task General
- Task MLC
- Task Planar
- Clinac 345
- Routine check
- Task CBCT
- Task Dosi
- Task General
- Task MLC
- Task Planar
- Clinac IX 34576
- TB STX 1023
- Dept Morning QA tests

Test Results Trend Analysis Patients Logoff

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
RAZOR^{DETECTOR}

High-performance diode detector


- ✓ Chip size: 0.95 x 0.95 mm; t = 0.5 mm
- ✓ Sensitive area: \varnothing 0.6 mm
- ✓ For Photon and Electron beams in RT

RAZOR^{CHAMBER}

Compact air ionization chamber


- ✓ Cavity volume: 0.01 ccm
- ✓ Central electrode material: graphite
- ✓ For Photon and Electron beams in RT

RAZOR^{NANOCHAMBER}

Smallest available ionization chamber

- ✓ Cavity volume: 0.003 ccm
- ✓ Central electrode material: graphite
- ✓ For Photon and Electron beams in RT

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